



GUIDE TO ISOLATION GOWN GUIDELINES

TO SELECT THE APPROPRIATE ISOLATION GOWN, IT'S HELPFUL TO UNDERSTAND THE STANDARDS AND TESTING REQUIREMENTS FOR AAMI - RATED GOWNS.

AAMI

ANSI/AAMI PB70:2012 provides standards for **liquid barrier performance**, with ratings for different levels of gown protection. The AAMI Liquid Barrier standard addresses fluid protection in the areas most likely to come in direct contact with potentially infectious material.

Because the AAMI standard requires that **isolation gowns** be fully covered in the back, open-back gowns do not meet the AAMI standard, and cannot carry an AAMI rating. Both the fabric and all seams and attachments must meet AAMI standards. In addition, all AAMI Levels must meet a 4% AQL (Acceptable Quality Limit).



MODERATE FLUID BARRIER
PROTECTION Use for the widest range of applications. Where moderate fluid protection is indicated.

Spray Impact Penetration ≤ 1.0g Hydrostatic Pressure ≥ 50cm



MINIMAL TO LOW FLUID BARRIER
PROTECTION for use when fluid
exposure is expected to be low.
Spray Impact Penetration ≤ 1.0g

Spray Impact Penetration ≤ 1.0g Hydrostatic Pressure ≥ 20cm



MINIMAL TO LOW FLUID BARRIER PROTECTION

Spray Impact Penetration ≤ 4.5g

FOR DETAILED GUIDELINES, GO TO www.aami.org

1 Association for the Advancement of Medical Instrumentation (AAMI), Liquid barrier performance and classification of protective apparel and drapes intended for use in health care facilities (ANSI/AAMI PB70:2012), May 2012, pp 6-7 (§4.2.1-4.2.3) https://www.aami.org/docs/default-source/products_store/standards/pb70_1206.pdf?sfvrsn=9d7b9544_2

FDA

As a medical device company, we are regulated by the US Food and Drug Administration (FDA). The manufacturing sites for HALYARD* AAMI Gowns are registered with the FDA, and the FDA reviews and audits our manufacturing and quality systems to ensure we maintain FDA regulatory standards.

The FDA categorizes medical devices into one of three classes – Class I, II, or III – based on their risks and the regulatory controls necessary to provide assurance of safety and effectiveness.² Isolation gowns are classified as Class I Medical Devices.

FOR MORE GO TO www.FDA.gov

GOWN TESTING

To make sure HALYARD* isolation gowns meet the standards set out by these organizations, they undergo a variety of standard tests, including:



FLUID AND BACTERIOPHAGE BARRIER

(ASTM 1671) Measures resistance of materials used in protective clothing to bloodborne pathogens using viral penetration at 2psig and ambient pressure.



IMPACT PENETRATION TEST (AATCC 42)

Measures resistance of fabrics to penetration of water by spray impact, as measured by weight gain of a blotter in grams.



RESISTANCE TO BLOOD PENETRATION

(ASTM F1670) Evaluates the resistance of a material to blood strikethrough, using synthetic blood with a surface tension and viscosity representative of blood and some body fluids.



HYDROSTATIC PRESSURE TEST (AATCC 127)

Tests for fluid resistance by measuring the force required for water to penetrate gown fabric and seams.

² U.S Food and Drug Administration (FDA), Premarket Notification Requirements Concerning Gowns Intended for Use in Health Care Settings, Guidance for Industry and Food and Drug Administration Staff, Doc. #1500025, pp 1-9 https://search.usa.gov/search?query=15000256affiliate=fda1.

HALYARD* AAMI RATED ISOLATION GOWNS

AAMI-RATED ISOLATION GOWN PORTFOLIO

HALYARD* Tri-Layer AAMI3 Isolation Gowns

| Code | Description | Neck Closure | Cuff | Size | Items/Case |
|-------|--------------------------------------|---------------|-----------------------|-----------|------------|
| 50697 | Over-The-Head Isolation Gown, Yellow | Over-the-Head | Thumb Hooks & Elastic | Universal | 100 |
| 50698 | Over-The-Head Isolation Gown, Yellow | Over-the-Head | Thumb Hooks & Elastic | X-Large | 100 |
| 54310 | Isolation Gown, Blue | Tape-Tab | Elastic | Large | 100 |
| 54311 | Isolation Gown, Blue | Tape-Tab | Elastic | X-Large | 100 |

HALYARD* Tri-Layer AAMI2 Isolation Gowns

| Code | Description | Neck Closure | Cuff | Size | Items/Case |
|-------|--------------------------------------|---------------|-----------------------|-----------|------------|
| 43008 | BASICS* Isolation Gown, Yellow | Over-the-Head | Thumb Hooks & Elastic | Large | 100 |
| 43009 | BASICS* Isolation Gown, Yellow | Over-the-Head | Thumb Hooks & Elastic | X-Large | 100 |
| 43146 | Over-the-Head Isolation Gown, Yellow | Over-the-Head | Thumb Hooks | Large | 100 |
| 43147 | Over-the-Head Isolation Gown, Yellow | Over-the-Head | Thumb Hooks | X-Large | 100 |
| 69979 | Isolation Gown, Yellow | Tape-Tab | Elastic | Large | 100 |
| 35460 | Isolation Gown, Yellow | Tape-Tab | Elastic | Large | 100 |
| 69988 | Isolation Gown, Yellow | Tape-Tab | Elastic | X-Large | 100 |
| 69986 | Isolation Gown, Yellow | Tape-Tab | Elastic | Universal | 100 |
| 69981 | Isolation Gown, Blue | Tape-Tab | Elastic | Universal | 100 |
| 69987 | Isolation Gown, Blue | Tape-Tab | Elastic | X-Large | 100 |
| 13962 | BASICS* Isolation Gown, Yellow | Tape-Tab | Elastic | Universal | 100 |
| 13961 | BASICS* Isolation Gown, Yellow | Tape-Tab | Elastic | X-Large | 100 |

HALYARD* Tri-Layer AAMI1 Isolation Gowns

| Code | Description | Neck Closure | Cuff | Size | Items/Case |
|-------|------------------------|--------------|---------|---------|------------|
| 54100 | Isolation Gown, Yellow | Tape-Tab | Elastic | Large | 100 |
| 54101 | Isolation Gown, Yellow | Tape-Tab | Elastic | X-Large | 100 |
| 54110 | Isolation Gown, Blue | Tape-Tab | Elastic | Large | 100 |

HALYARD* IMPERVIOUS **GOWNS**

HALYARD* PROCEDURE GOWNS TESTED FOR USE WITH 53 CHEMOTHERAPY DRUGS

These gowns meet the USP <800> guidelines for safe handling of hazardous drugs (HDs).[†]

Materials: Polypropylene, with polyethylene coating

Exposure Anticipated: Potential for heavy fluid or cytotoxic drug contact in front

| Code # | Description | Cuff | Size | Items/Case |
|--------|--------------------------------|-----------|----------|------------|
| 47346 | CHEMOPLUS*, Full Back, Blue | Knit Cuff | Large | 100 |
| 47347 | CHEMOPLUS*, Full Back, Blue | Knit Cuff | XX-Large | 100 |
| 69606 | Over-the-Head, Open Back, Blue | Knit Cuff | Large | 100 |
| 37284 | Over-the-Head, Open Back, Blue | Knit Cuff | XX-Large | 100 |

HALYARD* IMPERVIOUS COMFORT GOWNS

Materials: SMS Film Laminate

| Code # | Description | Cuff | Size | Items/Case |
|--------|--------------------------------|-------------|----------|------------|
| 69600 | Over-the-Head, Open Back, Blue | Knit Cuff | Large | 100 |
| 69601 | Over-the-Head, Open Back, Blue | Knit Cuff | XX-Large | 100 |
| 69602 | Over-the-Head, Open Back, Blue | Thumb Hooks | Large | 100 |
| 69603 | Over-the-Head, Open Back, Blue | Thumb Hooks | XX-Large | 100 |

HALYARD* IMPERVIOUS GOWNS

Materials: Poly-Coated Spunbond

| Code # | Description | Cuff | Size | Items/Case |
|--------|----------------------------------|-------------|---------|------------|
| 47003 | Over-the Head, Open Back, Yellow | Thumb Hooks | Large | 100 |
| 47004 | Over-the Head, Open Back, Yellow | Thumb Hooks | X-Large | 100 |

HALYARD* IMPERVIOUS GOWNS

Materials: Laminated Spunbond

| Code # | Description | Cuff | Size | Items/Case |
|--------|--------------------------------|-------------|---------|------------|
| 46946 | Over-the-Head, Open Back, Blue | Thumb Hooks | Large | 100 |
| 46947 | Over-the-Head, Open Back, Blue | Thumb Hooks | X-Large | 100 |

HALYARD* IMPERVIOUS OPEN BACK GOWN

Materials: Protective, Tri-Layer fabric with polyethylene coating

| Code # | Description | Cuff | Size | Items/Case |
|--------|--------------------------------|------------|-------|------------|
| 69455 | Over-the-Head, Open Back, Blue | Knit Cuffs | Large | 60 |

HALYARD* FILM GOWN IN DISPENSER BOX

Materials: Polyethylene film

| Code # | Description | Cuff | Size | Items/Case |
|--------|---|-------------|-----------|------------|
| 69490 | Over-the-Head, Perforated Open Back, Blue | Thumb Hooks | Universal | 75 |
| 69307 | Over-the-Head, Perforated Open Back, Blue | Thumb Hooks | Universal | 120 |
| 69316 | Over-the-Head, Perforated Open Back, Blue | Thumb Hooks | X-Large | 60 |
| 48414 | Over-the-Head, Perforated Open Back, Blue | Thumb Hooks | X-Large | 75 |

[†] Gowns should be changed in accordance with USP <800> guidelines (every 2–3 hours or immediately after a spill or splash). Tested for 53 Chemo drugs per test procedure ASTM F739-12.

HALYARD* ISOLATION GOWNS

HALYARD* TRI-LAYER OVER-THE-HEAD ISOLATION GOWNS

- Designed for quick and easy donning and doffing
- Made of SMS
- Open back panel provides ventilation
- Available with thumb hooks to hold sleeves in place under gloves

| Code # | Description | Cuff | Size | Items/Case |
|--------|---|--------------------------|-----------|------------|
| 32500 | Medium Weight, Tri-Layer, Open Back, Yellow | Thumb Hooks with Elastic | Universal | 100 |
| 32501 | Medium Weight, Tri-Layer, Open Back, Yellow | Thumb Hooks with Elastic | X-Large | 100 |
| 29549 | Medium Weight, Open Back, Yellow | Thumb Hooks | Universal | 100 |
| 29516 | Medium Weight, Open Back, Yellow | Thumb Hook | XX-Large | 100 |
| 32574 | Medium Weight, Open Back, Yellow, Folded 11" x 5.5" | Thumb Hooks | Universal | 100 |
| 32572 | Medium Weight, Open Back, Yellow, Folded 11" x 5.5" | Thumb Hooks | XX-Large | 100 |
| 43148 | BASICS*, Tri-Layer, Open Back, Yellow | Thumb Hooks | Universal | 100 |
| 43149 | BASICS*, Tri-Layer, Open Back, Yellow | Thumb Hooks | X-Large | 100 |
| 32502 | Light Weight, Tri-Layer, Open Back, Yellow | Thumb Hooks | Universal | 100 |
| 32503 | Light Weight, Tri-Layer, Open Back, Yellow | Thumb Hooks | X-Large | 100 |

HALYARD* FULL BACK PROCEDURE GOWNS

- Full back gowns with waist ties and hook-and-loop neckline closures
- Made of SMS fabric

| Code # | Description | Cuff | Size | Items/Case |
|--------|---|------|-----------|------------|
| 69025 | Tri-Layer fabric, Full Back, Waist Ties, Blue | Knit | Universal | 100 |
| 69028 | Tri-Layer fabric, Full Back, Waist Ties, Blue | Knit | X-Large | 100 |

LIGHTWEIGHT COVER GOWNS

- Made of single-layer spunbond fabric
- Choice of open or closed back styles

| Code # | Description | Cuff | Size | Items/Case |
|--------|---|-------------|-----------|------------|
| 69105 | Over-the-Head, Open Back w/ Tie Waist, Yellow | Thumb Hooks | Universal | 100 |
| 69106 | Over-the-Head, Open Back w/ Tie Waist, Yellow | Thumb Hooks | X-Large | 100 |
| 69124 | Traditional, Full Back, w/ Tie Waist and Neck, Yellow | Elastic | Universal | 100 |
| 69125 | Traditional, Full Back, w/ Tie Waist and Neck, Yellow | Elastic | X-Large | 100 |
| 69127 | Traditional, Full Back, w/ Tie Waist and Neck, Blue | Elastic | Universal | 100 |
| 69129 | Traditional, Full Back, w/ Tie Waist and Neck, Blue | Elastic | X-Large | 100 |

HALYARD* COVERALLS & **SCRUB SUITS**

HALYARD* EXTRA PROTECTION COVERALLS

Materials: Protective, Tri-Layer fabric

Features: Back gather, raglan sleeve, elastic cuff, open ankle

| Code # | Color | Size | Items/Case |
|--------|-------|-----------|------------|
| 10075 | White | Large | 24 |
| 10085 | White | X-Large | 24 |
| 10095 | White | XX-Large | 24 |
| 10096 | White | XXX-Large | 24 |

HALYARD* PROTECTIVE COVERALLS

Materials: Protective, Tri-Layer fabric

Features: Raglan sleeve, open cuffs and ankle

| Code # | Color | Size | Items/Case |
|--------|-------|----------|------------|
| 75631 | Blue | Large | 24 |
| 75641 | Blue | X-Large | 24 |
| 75651 | Blue | XX-Large | 24 |

HALYARD* SPUNBOND COVERALLS

Materials: Spunbond

Features: Set-in-sleeve, knit cuff, open ankle

| Code # | Color | Size | Items/Case |
|--------|-------|-----------|------------|
| 85631 | White | Large | 24 |
| 85641 | White | X-Large | 24 |
| 85651 | White | XX-Large | 24 |
| 85661 | White | XXX-Large | 24 |

HALYARD* SCRUB SUITS

Materials: Protective, Tri-Layer fabric

Features: Shirts – round neck, chest pocket; Pants – waist

ties, side and rear pockets, reinforced in-seam

Shirts

| Code # | Color | Size | Items/Case |
|--------|-------------|----------|------------|
| 69701 | Medium Blue | Medium | 48 |
| 69702 | Medium Blue | Large | 48 |
| 69703 | Medium Blue | X-Large | 48 |
| 69704 | Medium Blue | XX-Large | 48 |

Pants

| Code # | Color | Size | Items/Case |
|--------|-------------|----------|------------|
| 69711 | Medium Blue | Medium | 48 |
| 69712 | Medium Blue | Large | 48 |
| 69713 | Medium Blue | X-Large | 48 |
| 69714 | Medium Blue | XX-Large | 48 |

HALYARD* BASICS* SCRUB SUITS

Materials: Tri-Layer fabric

Features: Shirts – V-neck, 1 chest pocket and 2 hip

pockets; Pants – drawstring or elastic waist

Shirts

| Code # | Color | Size | Items/Case |
|--------|-----------|----------|------------|
| 47786 | Dark Blue | Medium | 48 |
| 47787 | Dark Blue | Large | 48 |
| 47788 | Dark Blue | X-Large | 48 |
| 47789 | Dark Blue | XX-Large | 48 |

Pants - Draw String Waist

| Code # | Color | Size | Items/Case |
|--------|-----------|----------|------------|
| 47782 | Dark Blue | Medium | 48 |
| 47783 | Dark Blue | Large | 48 |
| 47784 | Dark Blue | X-Large | 48 |
| 47785 | Dark Blue | XX-Large | 48 |

Pants - Elastic Waist

| Code # | Color | Size | Items/Case |
|--------|-----------|----------|------------|
| 47790 | Dark Blue | Medium | 48 |
| 47791 | Dark Blue | Large | 48 |
| 47792 | Dark Blue | X-Large | 48 |
| 47793 | Dark Blue | XX-Large | 48 |

HALYARD* LAB COATS AND LAB JACKETS

HALYARD* UNIVERSAL PRECAUTIONS LAB COATS

Materials: Protective Tri-Layer fabric

Features: Traditional collar, raglan sleeve, knit cuffs, 6 snaps, raglan sleeves, 3 reinforced pockets, back vent,

neck loop for hanging, size/name label

| Code # | Color | Size | Items/ Case |
|--------|-------|----------|----------------|
| 10040 | White | Small | 25 |
| 10041 | White | Medium | 25 |
| 10042 | White | Large | 25 |
| 10043 | White | X-Large | 25 |
| 10044 | White | XX-Large | 10 |
| 10045 | Blue | Small | 25 |
| 10046 | Blue | Medium | 25 |
| 10047 | Blue | Large | 25 |
| 10048 | Blue | X-Large | 25 |
| 10049 | Blue | XX-Large | 10 |

HALYARD* BASICS* PLUS LAB COATS

Materials: Protective, Tri-Layer fabric

Features: Knit collar, raglan sleeve, knit cuffs, 5 snaps, 3 reinforced pockets, neck loop for hanging, size/name label

| Code # | Color | Size | Items/Case |
|--------|-------|----------|------------|
| 10020 | White | Small | 25 |
| 10021 | White | Medium | 25 |
| 10022 | White | Large | 25 |
| 10023 | White | X-Large | 25 |
| 10024 | White | XX-Large | 25 |
| 10030 | Blue | Small | 25 |
| 10031 | Blue | Medium | 25 |
| 10032 | Blue | Large | 25 |
| 10033 | Blue | X-Large | 25 |
| 10034 | Blue | XX-Large | 25 |

HALYARD* BASICS* LAB COATS

Materials: Protective, Tri-Layer fabric

Features: Traditional collar, raglan sleeve, knit cuffs, 4

snaps, 2 reinforced pockets, size/name label

| Code # | Color | Size | Items/ Case |
|--------|-------|-----------|----------------|
| 10120 | White | Small | 25 |
| 10121 | White | Medium | 25 |
| 10122 | White | Large | 25 |
| 10123 | White | X-Large | 25 |
| 10124 | White | XX-Large | 25 |
| 10157 | White | XXX-Large | 25 |

HALYARD* SPUNBOND LAB COATS

Materials: Spunbond

Features: Traditional collar, set-in-sleeve, knit cuff,

5 snaps, 3 pockets

| Code # | Color | Size | Items/ Case |
|--------|-------|----------|----------------|
| 10231 | White | Small | 25 |
| 10232 | White | Medium | 25 |
| 10233 | White | Large | 25 |
| 10234 | White | X-Large | 25 |
| 10235 | White | XX-Large | 25 |

HALYARD* UNIVERSAL PRECAUTIONS LAB JACKETS

Materials: Protective, Tri-Layer fabric

Features: Knit collar, knit cuffs, 5 snaps, 3 reinforced pockets, neck loop for hanging, size/name label

| Code # | Color | Size | Items/ Case |
|--------|-------|---------|----------------|
| 10069 | White | Small | 25 |
| 10071 | White | Medium | 25 |
| 10072 | White | Large | 25 |
| 10073 | White | X-Large | 25 |

HALYARD* PROFESSIONAL JACKETS

Materials: Protective, Tri-Layer fabric

Features: Knit collar, knit cuffs, 5 snaps, 2 reinforced pockets, neck loop for hanging, size/name label

| Code # | Color | Size | Items/ Case |
|--------|-------|---------|----------------|
| 10067 | Blue | Medium | 24 |
| 10078 | Blue | Large | 24 |
| 10087 | Blue | X-Large | 24 |

HALYARD* MEDIUM WEIGHT SMS LAB JACKET

Materials: Protective, Tri-Layer SMS fabric

Features: Knit collar, set-in-sleeve, knit cuff, 4 snaps, 3

pockets

| Code # | Color | Size | Items/ Case |
|--------|-------|----------|----------------|
| 10334 | Blue | Small | 25 |
| 10335 | Blue | Medium | 25 |
| 10336 | Blue | Large | 25 |
| 10337 | Blue | X-Large | 25 |
| 10338 | Blue | XX-Large | 25 |

HALYARD* HEADWEAR

HALYARD* SPUNBOND BOUFFANT CAPS

Materials: Spunbond fabric

| Code # | Description | Size | Items/Case |
|--------|-----------------------------|-------------|----------------|
| 69801 | Spunbond Bouffant Cap, Blue | Large-24" | 100/bx, 5bx/cs |
| 69803 | Spunbond Bouffant Cap, Blue | X-Large-27" | 100/bx, 5bx/cs |
| 47322 | Spunbond Bouffant Cap, Red | Large-24" | 100/bx, 5bx/cs |
| 47323 | Spunbond Bouffant Cap, Red | X-Large-27" | 100/bx, 5bx/cs |

HALYARD* TRI-LAYER BOUFFANT CAPS

Materials: Tri-Layer fabric

| Code # | Description | Size | Items/Case |
|--------|--------------------------------------|-------------|----------------|
| 69083 | Tri-Layer fabric Bouffant Cap, White | Large-24" | 100/bx, 3bx/cs |
| 69086 | Tri-Layer fabric Bouffant Cap, White | X-Large-27" | 100/bx, 3bx/cs |
| 69088 | Tri-Layer fabric Bouffant Cap, Blue | Large-24" | 100/bx, 3bx/cs |

HALYARD* SURGICAL CAP

Materials: 4-Ply tissue with polyester reinforcements

| Code # | Description | Size | Items/Case |
|--------|--------------------|-----------|----------------|
| 69240 | Surgical Cap, Blue | Universal | 100/bx, 3bx/cs |

HALYARD* PROTECTIVE SURGICAL CAP

Materials: Tri-Layer fabric on sides, Spunbond fabric on top

| Code # | Description | Size | Items/Case |
|--------|-------------------------------|-----------|----------------|
| 69520 | Protective Surgical Cap, Blue | Universal | 100/bx, 3bx/cs |

HALYARD* COVER MAX* WITH CAPTURE TECHNOLOGY

Materials: Absorbent non-woven

Features: Expanded head coverage with option of deployable bouffant for added hair coverage

| Code # | Description | Size | Items/Case |
|--------|---|-----------|----------------|
| 47317 | Cover Max Surgical Cap with Deployable Bouffant | Universal | 100/bx, 3bx/cs |
| 47318 | Cover Max Surgical Cap | Universal | 100/bx, 3bx/cs |

HALYARD* SURGICAL HOOD

Materials: Tri-Layer fabric on sides, Spunbond fabric on top

| Code # | Description | Size | Items/Case |
|--------|--|-----------|----------------|
| 69110 | Protective Surgical Hood with Tie Neck, Blue | Universal | 100/bx, 3bx/cs |

HALYARD* FOOTWEAR

Shoe covers protect the patient from bacteria on the healthcare professional's shoes.

BASICS* SHOE COVERS

Materials: Spunbond fabric

Coverage: Shoe

| Code # | Description | Size | Items/Case |
|--------|--------------------|-----------|----------------|
| 69114 | BASICS* Shoe Cover | Universal | 100/bx, 3bx/cs |
| 69121 | BASICS* Shoe Cover | X-Large | 100/bx, 3bx/cs |

HALYARD* HEAVY-DUTY SHOE COVER

Materials: Tri-Layer fabric

Coverage: Shoe

| Code # | Description | Size | Items/Case |
|--------|-----------------------|-----------|----------------|
| 69551 | Heavy-Duty Shoe Cover | Universal | 100/bx, 3bx/cs |

X-TRA TRACTION* SHOE COVERS

Materials: Tri-Layer fabric; adhesive traction strips

Coverage: Shoe

| Code # | Description | Size | Items/Case |
|--------|----------------------------|-----------|----------------|
| 69252 | X-TRA TRACTION* Shoe Cover | Universal | 100/bx, 3bx/cs |
| 69254 | X-TRA TRACTION* Shoe Cover | X-Large | 80/bx, 3bx/cs |

ANKLE-GUARD* SHOE COVERS

Materials: Tri-Layer fabric; adhesive traction strips Coverage: Shoe and ankle, elastic at the ankle

| Code # | Description | Size | Items/Case |
|--------|-------------------------|-----------|---------------|
| 69253 | ANKLE-GUARD* Shoe Cover | Universal | 90/bx, 3bx/cs |
| 69353 | ANKLE-GUARD* Shoe Cover | X-Large | 90/bx, 3bx/cs |

HI-GUARD* REGULAR COVERAGE BOOTS

Materials: Tri-Layer fabric with film coating on foot area; adhesive traction strips

Features: Elastic at ankle and top of boot

Coverage: Shoe and lower leg

| Code # | Description | Size | Items/Case |
|--------|--------------------------------------|-----------|---------------|
| 69571 | HI-GUARD* Regular Full-Coverage Boot | Universal | 50/bx, 3bx/cs |
| 69671 | HI-GUARD* Regular Full-Coverage Boot | X-Large | 50/bx, 3bx/cs |

HI-GUARD* ULTRA COVERAGE BOOTS

Materials: Plastic film coating over Tri-Layer fabric; two 1" urethane foam traction strips

Features: Tie at ankle and elastic at top of boot

Coverage: Shoe and lower leg

| Code # | Description | Size | Items/Case |
|--------|-------------------------------|-----------|---------------|
| 69572 | HI-GUARD* Ultra Coverage Boot | Universal | 30/bx, 4bx/cs |
| 69672 | HI-GUARD* Ultra Coverage Boot | X-Large | 30/bx, 4bx/cs |

Caution: Footwear users should be aware that a risk of slipping exists. Hospital floor may be waxed or wet. Good technique and OSHA require that footwear be removed when contaminated or upon leaving the work area. Do not use nonconductive shoe covers in locations where flammable anesthetics are administered.

Note: Shoe dimensions can vary from street shoes to athletic shoes. When selecting a shoe cover size, please consider the outside dimensions of the shoe to ensure a proper fit.

For more information, please visit: www.halyardhealth.com

Call 1-844-HALYARD (1-844-425-9273) in the United States and Canada.



Customer-centric solutions include the tools, resources and support to help you provide the best in patient care.

KNOWLEDGE NETWORK* Accredited Education In-service Training and Technical Support Credentialed Sales Representatives Tools & Best Practices Clinical Research Commitment to Excellence Infection Prevention Initiatives Sustainability Leadership Heritage of Innovation

